

| | |
|---------------|---|
| Product Name: | LRP-WF-DDL-18 (-80) |
| Description: | <p>WiFi Digital Data Logger (DDL)</p> <p>This wireless digital data logger thermometer connects easily to WiFi network and provides continuous data transmission to the cloud-based interface. Designed for Ultra-Low (ULT) freezers and monitors temperatures below -80°C ensuring critical samples are not compromised.</p> <p>Features:</p> <ul style="list-style-type: none">• Complies with current CDC, VFC, and CFR Guidelines• Audible and visual alarms high/low temperature alarms• Includes 3 years Certification of Calibration• $^{\circ}\text{F}/^{\circ}\text{C}$ Switchable readout• Min/Max historical temperature point• 1 stainless steel detachable probe• Low battery indicator• User programmable logging interval• Stores up to 3 years of data (30 minute intervals)• Temperature data downloadable via WiFi network• Includes an FEP cable.• Operating Temperature Range: -90°C to 105°C (-130°F to 221°F)• Accuracy: $\pm 0.1^{\circ}\text{C}$ ($\pm 0.3^{\circ}\text{F}$) |