

Product Name:	LRP-WF-DDL-18 (-80)
Description:	<p data-bbox="441 268 1432 304">WiFi Digital Data Logger (DDL)</p> <p data-bbox="441 340 1432 562">This wireless digital data logger thermometer connects easily to WiFi network (2.4GHz Radio Frequency) and provides continuous data transmission to the cloud-based interface. Designed for Ultra-Low (ULT) freezers and monitors temperatures below -80°C ensuring critical samples are not compromised. Features:</p> <ul data-bbox="483 598 1432 1108" style="list-style-type: none"> • Complies with current CDC, VFC, and CFR Guidelines • Audible and visual alarms high/low temperature alarms • Includes 2 years Certification of Calibration • $^{\circ}\text{F}/^{\circ}\text{C}$ Switchable readout • Min/Max historical temperature point • 1 stainless steel detachable probe • Low battery indicator • User programmable logging interval • Stores up to 3 years of data (30 minute intervals) • Temperature data downloadable via WiFi network • Includes an FEP cable. • Operating Temperature Range: -90°C to 105°C (-130°F to 221°F) • Accuracy: $\pm 0.1^{\circ}\text{C}$ ($\pm 0.3^{\circ}\text{F}$)