

Product Name:	CRT-LRP-HC-S49S
Description:	<p>Upright Controlled Room Temperature Cabinet, Solid Door  49 Cu. Ft. capacity  Two swing solid doors, self-closing, right hinged, not reversible  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  Adjustable operating temperature range: 68°F–77°F (20°C – 25°C)  Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Keyed door locks  Magnetic door gasket for positive seal  2 probes; 1 in air and 1 in sample bottle  Eight adjustable shelves  Vaccine storage power cord warning label  Forced draft air circulation  Heavy duty door handles  LED interior lights are safety shielded and switch controlled  Swivel Casters (3")  probe access port (3/8")  White powder-coated steel exterior and interior  Hydrocarbon, natural refrigerant (R290)  One year parts and labor warranty, plus an additional four year compressor parts warranty  Exterior dimensions: 54" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters  Cycle defrost  115V, 60 Hz, 4.5 Amps, 1/4 HP  ENERGY STAR CERTIFIED  ETL listed  Shipping weight: 446 lbs.</p>