

1

TempLog Plus

2

Plus Series

3

Basic Series

LRP
LAB RESEARCH
PRODUCTS



- Touch screen display
- Digital microprocessor
- Digital temperature display
- Audible and visual alarms
- Digital 24-hour temp. chart
- Digital data logging
- USB port for data transfer
- 4-20mA output
- Remote alarm contacts
- Door ajar alarm
- Low battery alarm
- Operation indicator
- Temp. scale bar display
- Sample temperature display

- Bi-Color LED display
- Digital microprocessor
- Digital temperature display
- Audible and visual alarms
- Resettable Min/Max
- Remote alarm contacts
- Power failure alarm*
- Sensor error alarm
- Alarm mute w/ring-back*
- Sample temperature display

**Not available on all models*
See specs for additional information*

- Digital temperature display*
- Digital microprocessor*
- Audible and visual alarms*
- Resettable Min/Max*
- Remote alarm contacts*

**Not available on all models*
See specs for additional information*



The TempLog Plus Series refrigerators and freezers are purpose built for scientific and medical applications requiring superior temperature management and state of the art product security and data logging.

The Plus Series refrigerators and freezers are also purpose built for critical cold storage applications, but where built-in data logging functionality may not be a necessity.

The Basic Series refrigerators and freezers can help meet the demand of scientific and medical applications at budget-easing price points or satisfy any general purpose cold storage needs.