

|               |   |
|---------------|---|
| Product Name: | LRP-HC-UCBI-0404G-ADA   |
| Description:  | <p>Plus Series Undercounter Refrigerator Built-In<br/> 4.6 Cu. Ft. Capacity<br/> One glass swing door, right hinged*<br/> Microprocessor temperature controller<br/> Digital temperature display<br/> Audible and visual high and low temperature alarms<br/> Remote alarm contacts<br/> 2 probes; 1 in air and 1 in glycol bottle<br/> Adjustable operating temperature range: 1°C – 10°C<br/> Three shelves (two adjustable/one fixed) with guard rail on back to prevent product from pushing against back wall<br/> Back wall evaporator cover<br/> Forced draft circulation<br/> Keyed door lock<br/> Magnetic door gasket for positive seal<br/> Heavy duty door handle<br/> Leveling legs<br/> LED interior light<br/> Probe access port (3/8")<br/> Powder coat exterior<br/> High density urethane foam cabinet and door insulation<br/> HFC-free refrigerant (R600a), foam insulation and packaging<br/> Two year parts and labor warranty, plus an additional three year compressor parts only warranty<br/> Pyxis®, Omnicell® and AcuDose RX® compatible<br/> Exterior dimensions: 23 ¾ "W x 25 ½" D x 31 15/16 "H<br/> Cycle defrost<br/> 115V, 60 Hz, 1.74 Amps, 1/8 HP<br/> ENERGY STAR CERTIFIED<br/> Energy Consumption 1.15 kWh/day<br/> ADA compliant<br/> UL/C-UL listed<br/> Shipping weight: 150 lbs<br/> *This product is available left hinged, product number LRP-HC-UCBI-0404G-ADA-LH</p> |