

Product Name:	LRP-HC-UCFS-0504
Description:	<p>Plus Series Undercounter Refrigerator Freestanding  5.2 Cu. Ft. Capacity  One solid swing door, right hinged, not reversible  Microprocessor temperature controller  Adjustable operating temperature range: 1°C – 10°C  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  2 probes; 1 in air and 1 in glycol bottle  Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall  Bottom footed shelf to allow air circulation  Back wall evaporator cover  Forced draft circulation  Keyed door lock  Magnetic door gasket for positive seal  Leveling legs  LED interior light  Switch activated light  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet insulation  HFC-free refrigerant (R600a), foam insulation and packaging Pyxis<sup>®</sup>, Omnicell<sup>®</sup> and AcuDose RX<sup>®</sup> compatible  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Exterior dimensions: 23 3/4"W x 24" D x 32 1/8 "H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Cycle defrost  115V, 60 Hz, 1.3 Amps, 1/6 HP  ENERGY STAR CERTIFIED  UL/C-UL listed  Shipping weight: 129 lbs</p>