

Product Name:	LRP-HC-RFC20A
Description:	<p>Basic Series Refrigerator & Freezer Combination 20 cu. ft. capacity (15.35 ref / 5.15 freezer) Two exterior solid swing doors, right hinged Mechanical thermostat Refrigerator temperature range : 0°C to 10°C Freezer chamber: -15°C to -10°C with two auto defrost cycles warming to 0°C per day Three shelves, refrigerator (one fixed, two adjustable); one fixed shelf, freezer Magnetic door gasket for positive seal Heavy duty door handles Leveling legs Interior refrigerator light Two probe access ports (3/8"); one, refrigerator/one, freezer White textured cabinet and door High density urethane foam cabinet and door insulation Hydrocarbon, natural refrigerant (R600a) Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 30"W x 34 1/8"D x 69"H Auto defrost (Freezer); Cycle Defrost (Refrigerator) 115V, 60 Hz, 3 Amps Shipping weight: 231 lbs. DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assist in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.</p>