

Product Name:	RFC-30M
Description:	<p>Refrigerator & Freezer Combination Basic Series 30 Cu. Ft. capacity (25 Ref./5 Freezer) One exterior solid swing door and one interior solid swing door, right hinged Mechanical thermostat Refrigerator temperature range: 0°C to 10°C Freezer temperature range: -15°C to -20°C 3 adjustable shelves, refrigerator Keyed door lock Magnetic door gasket for positive seal Stainless steel door handle Leveling Legs Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer) Powder coat paint exterior High density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging Accu-test certified One year parts and labor warranty including compressor Exterior dimensions: 35 ¾ "W x 35 ½ "D x 72 ½ "H Manual defrost 115V, 60 Hz, 6 Amps, 1/3 HP Shipping weight: 383 lbs.</p>