

Product Name:	LRP-HC-RFC7S
Description:	<p>Plus Series Refrigerator &amp; Freezer Combination</p> <p>5.2 cu. ft. capacity (Refrigerator) / 1.7 cu. ft. capacity (Freezer)</p> <p>Two exterior swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible</p> <p>Microprocessor temperature controller, refrigerator and freezer</p> <p>Digital temperature display, refrigerator and freezer</p> <p>Audible and visual high and low temperature alarms, refrigerator and freezer</p> <p>2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer</p> <p>Remote alarm contacts, refrigerator and freezer</p> <p>Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer:-15°C – -25°C</p> <p>3 adjustable shelves (refrigerator) / 2 adjustable shelves (freezer), shelves have back rails to prevent product from being pushed against back wall</p> <p>Footed bottom shelf refrigerator, and freezer</p> <p>Keyed door locks</p> <p>Magnetic door gasket for positive seal</p> <p>Integrated door handles</p> <p>Leveling legs</p> <p>Stacking kit</p> <p>Two probe access ports (3/8"); one, refrigerator/one, freezer White powder coated exterior</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Two year parts and labor warranty, plus an additional three year compressor parts warranty</p> <p>Exterior dimensions: 23 3/4"W x 24"D x 53 1/8"H</p> <p>Manual defrost freezer / cycle defrost refrigerator</p> <p>115V, 60 Hz, 1.3 Amps refrigerator, 1 Amp freezer, 1/6 HP refrigerator, 1/8 HP freezer</p> <p>UL/C-UL listed</p> <p>Shipping weight: 215 lbs.</p> <p>*Please note this unit requires 2 outlets.</p>