

Product Name:	LRP-HC-RFC7SA
Description:	<p>Plus Series Refrigerator & Auto Defrost Freezer Combination 5.2 cu. ft. (Refrigerator) / 1.3 cu. ft. capacity (Freezer) Two exterior swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible Microprocessor temperature controller, refrigerator and freezer Digital temperature display, refrigerator and freezer Audible and visual high and low temperature alarms, refrigerator and freezer 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer Remote alarm contacts, refrigerator and freezer Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer:-15°C – -25°C 3 adjustable shelves (refrigerator)/2 adjustable shelves (freezer), shelves have back rails to prevent product from being pushed against back wall Footed bottom shelf refrigerator, and freezer Keyed door locks Magnetic door gasket for positive seal Integrated door handles Leveling legs Stacking kit Two probe access ports (3/8"); one, refrigerator/one, freezer White powder coated exterior HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 23 3/4"W x 24"D x 53 1/8"H Auto defrost freezer / cycle defrost refrigerator 115V, 60 Hz, 1.3 Amps refrigerator, 1 Amps freezer, 1/6 HP refrigerator, 1/8 HP freezer UL/C-UL listed Shipping weight: 215 lbs. *Please note this unit requires 2 outlets.</p>