

Product Name:	LRP-HCRFC9SS-LH
Description:	<p>Plus Series Refrigerator &amp; Freezer Combination-LEFT HINGED  9 Cu. Ft. Capacity (4.5 Ref./4.2 Freezer)  Two stainless steel exterior swing doors, left hinged  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 1°C – 10°C;  Freezer:-15°C – -25°C  Audible and visual high and low temperature alarms, refrigerator  and freezer  2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and  Freezer  Remote alarm contacts, refrigerator and freezer  2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with  guard rails on back in Refrigerator to prevent product from being  pushed against back wall  Back wall evaporator cover, Refrigerator  Forced draft circulation Refrigerator  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handles  Leveling legs  Stacking kit  LED interior lights  Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer)  Stainless steel exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Two year parts and labor warranty, plus an additional three year  compressor parts warranty  Exterior dimensions: 23 ¾"W x 24 ½"D x 66 ¾"H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP  Refrigerator and Freezer  UL/C-UL listed  Shipping weight: 266 lbs.  *Please note this unit requires 2 outlets.</p>