

Product Name:	LRP-FFB-30
Description:	<p>Basic Series Flammable Storage Freezer</p> <p>30 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Mechanical thermostat</p> <p>Adjustable operating temperature range: -15°C – -25°C</p> <p>Four fixed shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Stainless steel door handle</p> <p>Leveling legs</p> <p>Probe access port (3/8")</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>CFC-free refrigerant (R404a), foam insulation and packaging</p> <p>One year parts and labor warranty</p> <p>Exterior dimensions: 35 ¾ "W x 35 ½ "D x 72 ½ "H</p> <p>Manual defrost</p> <p>115V, 60Hz, 6 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Units meet NFPA guidelines 45 and 70</p> <p>Shipping weight: 378 lbs.</p>