

Product Name:	LRP-FRP-04
Description:	<p>Basic Series Undercounter Flammable Storage Refrigerator</p> <p>Freestanding 5 Cu. Ft. capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Microprocessor temperature controller</p> <p>2 probes; 1 in air and 1 in sample bottle</p> <p>Standard digital temperature display</p> <p>Audible and visual high and low temperature alarms</p> <p>Remote alarm contacts</p> <p>Adjustable operating temperature range: 1°-10°C</p> <p>Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Integrated door handle</p> <p>One probe access port (3/8")</p> <p>White cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Two year parts and labor warranty</p> <p><i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i></p> <p>Exterior dimensions: 23 5/8" W x 24 1/4" D x 32 3/16" H</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 1.2 Amps, 1/9 HP</p> <p>UL/C-UL listed</p> <p>Units meet NFPA guidelines 45 and 70</p> <p>Shipping weight: 114 lbs.</p>