

Product Name:	LRP-HC-EFP-17
Description:	<p>Basic Series Hazardous Location (Explosion Proof) Freezer 17 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Microprocessor Temperature Controller Adjustable operating temperature range: -15°C – -25°C Six fixed shelves Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle Leveling legs White textured cabinet and door High density urethane foam cabinet and door insulation Hydrocarbon, natural refrigerant (R290)</p> <p>Exterior dimensions: 30 3/4"W x 28 5/8"D x 67 3/4"H Manual defrost Not equipped with plug; must be hardwired in conduit to 120 VAC power source 115V, 60Hz, 5.5 Amps, 1/2 HP ETL listed Two year parts and labor warranty</p> <p>Meet NFPA and OSHA guidelines for hazardous locations, Class 1, Division II, Group C & D. These units are engineered for use in environments where volatile / explosive conditions could potentially exist. For storage of flammable or volatile materials within typical laboratory or non-volatile environments, see Flammable Material Storage units. Shipping weight: 220 lbs.</p>