

Product Name:	LRP-HC-LB-69
Description:	<p>Basic Series Glass Door Laboratory Refrigerator  69 Cu. Ft. Capacity  Three sliding glass doors, self-closing (energy efficient, double pane vacuum insulated)  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature: 1°C to 10°C  Twelve adjustable shelves (adjustable in ½" increments)  Forced draft air circulation  Keyed door lock  Magnetic door gasket for positive seal  Integrated door handles  Swivel Casters (3 ½")  LED interior lights are safety shielded and switch controlled  Probe access port (¾")  Powder coated steel white interior and exterior  High density urethane foam cabinet insulation  Hydrocarbon, natural refrigerant (R290)  One year parts and labor warranty, plus an additional four year compressor parts warranty  Exterior dimensions: 81" W x 33 ¼" D x 81 ¾" H with casters; 78 ¼" H without casters  Cycle defrost  115V, 60 Hz, 7.5 Amps, 1/2 HP  ETL listed  Shipping weight: 738 lbs.</p>