

Product Name:	LRP-HC-RFP-17
Description:	<p>Basic Series Natural Refrigerant Laboratory Refrigerator  17 Cu. Ft. capacity  One solid swing door, right hinged, not reversible  Microprocessor Temperature Controller  Adjustable operating temperature range: 1°-10°C  Digital Temperature Display  Sample simulated digital temperature display; .1° Resolution  Audible and visual high/low temperature alarms  Remote alarm contacts  5 adjustable shelves  Keyed door lock  Magnetic door gasket for positive seal  Heavy duty door handle  Leveling legs  One probe access port (3/4")  White textured cabinet and door  High density urethane foam cabinet and door insulation  Natural, Hydrocarbon Refrigerant (R600a)  One year parts and labor warranty, plus an additional four-year  compressor parts warranty  Exterior dimensions: 30 3/4" W x 28 5/8"D x 67 3/4"H  Cycle defrost  115V, 60 Hz, 2.5 Amps, 1/4 HP  ENERGY STAR CERTIFIED  ETL listed  Shipping weight: 240 lbs.</p>