

Product Name:	LRP-HC-RFP-17
Description:	<p>Basic Series Natural Refrigerant Laboratory Refrigerator 17 Cu. Ft. capacity One solid swing door, right hinged, not reversible Microprocessor Temperature Controller Adjustable operating temperature range: 1°-10°C Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms Sample simulated digital temperature display; .1° Resolution Adjustable temperature display in Celsius or Fahrenheit Audible and visual high/low temperature alarms Sensor error alarm Power failure alarm Resettable 1-point Min/Max History Remote alarm contacts Alarm validation (audible, visual, and remote alarm) Alarm mute with ring back 5 adjustable shelves Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle Leveling legs One probe access port (3/4") White textured cabinet and door High density urethane foam cabinet and door insulation Natural, Hydrocarbon Refrigerant (R600a) One year parts and labor warranty, plus an additional four-year compressor parts warranty Exterior dimensions: 30 3/4" W x 28 5/8"D x 67 3/4"H Cycle defrost 115V, 60 Hz, 3 Amps, 1/4 HP ETL listed Shipping weight: 240 lbs.</p>