

Product Name:	LRP-HC-MFP-17
Description:	<p>Basic Series Natural Refrigerant Manual Defrost Laboratory Freezer</p> <p>17 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged</p> <p>Microprocessor temperature controller</p> <p>Digital temperature display</p> <p>Audible and visual high and low temperature alarms Remote alarm contacts</p> <p>2 probes; 1 in air and 1 in glycol bottle</p> <p>Adjustable operating temperature range: -15°C – -25°C</p> <p>Six fixed shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Leveling legs</p> <p>Probe access port (3/4")</p> <p>Powder coat exterior</p> <p>High density urethane foam cabinet and door insulation</p> <p>Natural, Hydrocarbon Refrigerant (R290)</p> <p>Two year parts and labor warranty, plus additional three year compressor parts warranty</p> <p>Exterior dimensions: 30 3/4" W x 28 5/8" D x 67 3/4"H Note: This unit must have 4" clearance on sides and back for adequate ventilation</p> <p>Manual defrost</p> <p>115V, 60 Hz, 5.5 Amps, ½ HP</p> <p>ETL listed</p> <p>Shipping weight: 220 lbs.</p>