

Product Name:	LRP-MFB-20-C
Description:	<p>Basic Series Manual Defrost Laboratory Chest Freezer 20 Cu. Ft. capacity One solid door Mechanical thermostat Adjustable operating temperature range: -15°C – -25°C LED Interior light Two removable baskets Swivel casters (2 ¾") Keyed door lock Magnetic door gasket for positive seal White cabinet and lid One probe access port (3/8") High density urethane foam cabinet and door insulation CFC-free refrigerant (R134a), foam insulation and packaging One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 74 3/8"W x 29½"D x 34"H with casters; 31 ¼"H without casters (Height with hinged door open 90° is 61 ¾") Manual defrost 115V, 60 Hz, 5 Amps, 1/4 HP Shipping weight: 241 lbs.</p>