

Product Name:	LRP-MFP-15-C
Description:	<p>Plus Series Manual Defrost Laboratory Chest Freezer  15 Cu. Ft. capacity  One solid door  Microprocessor temperature controller:  Adjustable operating temperature range: -15°C – -25°C  Standard Digital Temperature Display and Alarm Module  Sample simulated digital temperature display: .1C resolution,  convertible to display temperatures in Celsius or Fahrenheit  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Resettable 1-point Min/Max History  Remote alarm contacts  LED Interior light  Two removable baskets  Swivel casters (2 ¾")  Keyed door lock  Magnetic door gasket for positive seal  White cabinet and lid  One probe access port (3/8")  High density urethane foam cabinet and door insulation  CFC-free refrigerant (R134a), foam insulation and packaging  Two year parts and labor warranty; plus additional three year  compressor parts warranty  Exterior dimensions: 56" W x 29 ½"D x 35 1/8"H with casters; 32  3/8" without casters  (Height with hinged door open 90° is 62 ½")  Manual defrost  115V, 60 Hz, 5 Amps, 1/4 HP  Shipping weight: 197 lbs.</p>