

Product Name:	LRP-MFP-20-C
Description:	<p>Plus Series Manual Defrost Laboratory Chest Freezer  20 Cu. Ft. capacity  One solid door  Microprocessor temperature controller:  Adjustable operating temperature range: -15°C – -25°C  LED Interior light  Two removable baskets  Swivel casters (2 ¾")  Standard Temperature Display &amp; Alarm Module:  Sample simulated digital temperature display: .1C resolution,  convertible to display temperatures in Celsius or Fahrenheit  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Resettable 1-point Min/Max History  Remote alarm contacts  Keyed door lock  Magnetic door gasket for positive seal  White cabinet and lid  One probe access port (3/8")  High density urethane foam cabinet and door insulation  CFC-free refrigerant (R134a), foam insulation and packaging  One year parts and labor warranty, plus an additional four year  compressor parts warranty  Exterior dimensions: 74 3/8"W x 29½"D x 34"H with casters; 31 ¼"H  without casters (Height with hinged door open 90° is 61 ¾")  Manual defrost  115V, 60 Hz, 5 Amps, 1/4 HP  UL/C-UL listed  Shipping weight: 241 lbs.</p>