

Product Name:	LRP-HCPTP-47
Description:	<p>Plus Series Pass Thru Laboratory Refrigerator 47 Cu. Ft. capacity Two sliding glass doors on front, two solid swing doors on back Microprocessor Temperature Controller Adjustable operating temperature range: 1°-10°C Digital Temperature Display and Alarm Module Sample simulated digital temperature display; .1° Resolution Adjustable temperature display in Celsius or Fahrenheit Audible and visual high/low temperature alarms Sensor error alarm Power failure alarm Resettable 1-point Min/Max History Remote alarm contacts Alarm validation (audible, visual, and remote alarm) Alarm mute with ring back Eight adjustable shelves Forced draft air circulation Magnetic door gasket for positive seal Heavy duty door handle Casters (adds 4" to height) One probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R290), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 54 1/8"W x 32 3/8"D x 82 5/8" with casters;78 5/8 "H without casters Cycle defrost 115V, 60 Hz, 8.8 Amps, ½ HP, BTU/HR 7315 UL/C-UL listed Shipping weight: 575 lbs.</p>