

Product Name:	PH-LRP-HC-UCBI-0420-ADA
Description:	<p>Plus Series Pharmacy/Vaccine Undercounter Freezer Built-In 4.2 Cu. Ft. Capacity One solid swing door, right hinged* Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in sample bottle Adjustable operating temperature range: 5°F – -13°F (-15°C – -25°C) Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory, F/C switchable, field installable, and visual & audible temp alarm Two fixed shelves Vaccine storage power cord warning label Keyed door lock Magnetic door gasket for positive seal Recessed door handle Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 23 ¾"W x 24 ½"D x 31 15/16"H Manual defrost 115V, 60 Hz, 1.5 Amps, 1/8 HP ADA compliant UL/C-UL listed ENERGY STAR CERTIFIED Shipping weight: 133 lbs. *This product is available left hinged, product number PH-LRP-HC-UCBI-0420-ADA-LH</p>