

Product Name:	LRP-HC-SSP-23
Description:	<p>Plus Series Stainless Steel Laboratory Refrigerator 23 Cu. Ft. Capacity One stainless steel swing door, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Three adjustable shelves Forced draft air circulation Keyed door lock Magnetic door gasket for positive seal Integrated door handle Swivel casters (adds 6 ½" to height) Interior light Probe access port (3/8") Stainless steel interior floor Stainless steel exterior High density urethane foam cabinet insulation Hydrocarbon, natural refrigerants (R290) Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 26 ¾"W x 34 ⅞"D x 83 ⅞"H with casters; 77 ⅜"H without casters Cycle defrost 115V, 60 Hz, 3 Amps, 1/5 HP ETL listed Shipping weight: 295 lbs.</p>