

Product Name:	LRP-HC-UCFS-0420W
Description:	<p>Basic Series Undercounter Freezer Freestanding 4 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat Adjustable operating temperature range: -15°C – -25°C 2 fixed shelves Keyed door lock Magnetic door gasket for positive seal Recessed door handle Leveling legs Probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging One year parts and labor warranty Exterior dimensions: 21 1/2"W x 22 1/4"D x 32 3/4"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost 115V, 60 Hz, .90 Amps, 1/7 HP UL/C-UL listed Shipping weight: 94 lbs.</p>