

Product Name:	LRP-HC-UCBI-0420-LH
Description:	<p>Plus Series Undercounter Built in Freezer-LEFT HINGED 4.2 Cu. Ft. Capacity One solid swing door, left hinged Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in glycol bottle Adjustable operating temperature range: -15°C – -25°C Two fixed shelves Keyed door lock Recessed door handle Magnetic door gasket for positive seal Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Exterior dimensions: 23 3/4"W x 24 1/2"D x 33 3/8 "H Manual defrost 115V, 60 Hz, 1.5 Amps 1/8 HP ENERGY STAR CERTIFIED Energy Consumption 0.97 kWh/day UL/C-UL listed Shipping weight: 132 lbs.</p>