

Product Name:	PH-LRP-HCRFC9SS
Description:	<p>Plus Series Pharmacy/Vaccine Refrigerator & Freezer Combination 9 Cu. Ft. Capacity (4.5 Ref./4.2 Freezer) Two stainless steel exterior swing doors, right hinged* Microprocessor temperature controller, refrigerator and freezer Digital temperature display, refrigerator and freezer Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C) Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable Audible and visual high and low temperature alarms, refrigerator and freezer Remote alarm contacts, refrigerator and freezer 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer 2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall Back wall evaporator cover, Refrigerator Forced draft circulation Refrigerator Vaccine storage power cord warning label Keyed door locks Magnetic door gasket for positive seal Recessed door handles Leveling legs Stacking kit LED interior lights Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer) Stainless steel exterior HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 23 3/4"W x 24 1/2"D x 66 3/4"H Manual defrost freezer / cycle defrost refrigerator 115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer UL/C-UL listed Shipping weight: 266 lbs. *Please note this unit requires 2 outlets. *This product is available left hinged, product number PH-LRP-HC-RFC9SS-LH</p>