

|               |  |
|---------------|--|
| Product Name: | PH-LRP-HC-S12S   |
| Description:  | <p>Basic Series Pharmacy/Vaccine Solid Door Refrigerator<br/> 12 Cu. Ft. Capacity<br/> One swing solid door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible<br/> Microprocessor temperature controller<br/> Digital temperature display<br/> Audible and visual high and low temperature alarms<br/> Remote alarm contacts<br/> 2 probes; 1 in air and 1 in sample bottle<br/> Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)<br/> Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual &amp; audible temp alarm<br/> Three adjustable shelves (adjustable in ½” increments)<br/> Vaccine storage power cord warning label<br/> Forced draft air circulation<br/> Keyed door lock<br/> Magnetic door gasket for positive seal<br/> Heavy duty door handle<br/> Swivel Casters (3 ½”)<br/> LED interior lights are safety shielded and switch controlled<br/> Probe access port (3/8”)<br/> Powder coated steel white interior and exterior<br/> High density urethane foam cabinet insulation<br/> Hydrocarbon, natural refrigerants (R290)<br/> Pharmacy refrigerator/freezer toolkit and temperature logs<br/> One year parts and labor warranty, plus an additional four year compressor parts warranty<br/> Exterior dimensions: 25” W x 29 ¾” D x 65 ¾” H with casters; 62 ¼” H without casters<br/> Cycle defrost<br/> 115V, 60 Hz, 3 Amps, 1/5 HP<br/> ETL listed<br/> Shipping weight: 227 lbs.</p> |