

Product Name:	PH-LRP-HC-UCBI-0204
Description:	<p>Plus Series Pharmacy/Vaccine Undercounter Refrigerator Built-In 2.5 Cu. Ft. Capacity One solid swing door, right hinged Microprocessor temperature controller Digital temperature display Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in glycol bottle Remote alarm contacts Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall Back wall evaporator cover Vaccine storage power cord warning label Forced draft circulation Self-closing door Keyed door lock Magnetic door gasket for positive seal Leveling legs LED light Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 17 3/4"W x 22 3/4"D x 30 3/4"H Energy Star Certified for Pharmacy Refrigerators and Freezers Cycle defrost 115V, 60 Hz, .90 Amps, 1/8 HP Energy Consumption 0.58 kWh/day ADA compliant UL/C-UL listed Shipping weight: 85 lbs.</p>