

Product Name:	LRP-HC-SSP-72G
Description:	<p>Plus Series Stainless Steel Glass Door Refrigerator 72 Cu. Ft. Capacity Three swing glass doors, self-closing, not reversible 12 ea. Adjustable Shelves Digital display microprocessor temperature controller Adjustable operating temperature: 1°C to 10°C (34°F to 50°F) Audible and visual Hi/Lo temperature alarm 2 probes; 1 in air and 1 in sample bottle Sensor error alarm Continuous product (simulated) temperature display Remote alarm contacts Forced draft air circulation Keyed door locks Magnetic door gaskets for positive seal Heavy duty door handles LED interior lights are safety shielded and switch controlled Probe access port (3/4") Stainless steel exterior White-powder coated aluminum interior with 304 stainless steel floor Easily cleaned smooth interior Hydrocarbon natural refrigerant (R290) 2 year parts and labor warranty, plus an additional 3 year compressor parts warranty 6 swivel casters (5 3/4") Exterior dimensions: 81" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters Cycle defrost 115V, 60 Hz, 7.5 Amps, 1/2 HP ETL listed Shipping weight: 786 lbs.</p>