

Product Name:	CRT-LRP-HC-S36S
Description:	Upright Controlled Room Temperature Cabinet, Solid Door 36 Cu. Ft. capacity Two swing solid doors, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 68°F–77°F (20°C – 25°C) Keyed door locks Magnetic door gasket for positive seal 2 probes; 1 in air and 1 in sample bottle Calibrated using a NIST traceable device, certificate included Eight adjustable shelves Vaccine storage power cord warning label Forced draft air circulation Heavy duty door handles LED interior lights are safety shielded and switch controlled Swivel Casters (3") Probe access port ($\frac{3}{8}$") White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 39 $\frac{5}{8}$" W x 34 $\frac{3}{4}$" D x 81 $\frac{3}{4}$" H with casters; 78 $\frac{1}{4}$" H without casters Cycle defrost 115V, 60 Hz, 4.5 Amps, 1/4 HP ETL listed Shipping weight: 360 lbs.