

Product Name:	PH-LRP-HC-RFC1020G
Description:	Plus Series Pharmacy Refrigerator & Freezer Combination (-20°C) Freestanding 5.2 cu. ft. capacity (Refrigerator) / 4 cu. ft. capacity (Freezer) One glass swing door (fridge), one solid swing door (freezer) Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range - Freezer: 5°F – -13°F (-15°C – -25°C); Refrigerator: 36°F to 46°F (2°C to 8°C) Calibrated using a NIST traceable device, certificate included 2 fixed shelves (freezer); 3 adjustable shelves (fridge) Bottom footed shelf to allow air circulation Back wall evaporator cover Keyed door locks Magnetic door gasket for positive seal Recessed door handles Stacking kit Leveling legs LED interior light Probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Exterior dimensions: 23 3/4" W x 24"D x 65"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost (freezer); Cycle defrost (refrigerator) 115V, 60 Hz 1.4 Amps, 1/8 HP (freezer); 1.3 Amps, 1/6 HP (refrigerator) UL/C-UL listed Shipping weight: 266 lbs.