

| | |
|---------------|---|
| Product Name: | PH-LRP-HC-RFC9G-LH |
| Description: | <p>Plus Series Pharmacy/Vaccine Refrigerator & Freezer Combination-LEFT HINGED</p> <p>9 Cu. Ft. Capacity (4.6 Ref./4.2 Freezer)</p> <p>Two exterior swing doors (glass door refrigerator/solid door freezer), left hinged</p> <p>Microprocessor temperature controller, refrigerator and freezer</p> <p>Digital temperature display, refrigerator and freezer</p> <p>Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)</p> <p>Calibrated using a NIST traceable device, certificate included</p> <p>Audible and visual high and low temperature alarms, refrigerator and freezer</p> <p>Remote alarm contacts, refrigerator and freezer</p> <p>2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer</p> <p>2 adjustable/1 fixed shelves refrigerator with guard rail on back to prevent product from being pushed against back wall/2 fixed shelves freezer</p> <p>Back wall evaporator cover, Refrigerator</p> <p>Forced draft circulation, Refrigerator</p> <p>Vaccine storage power cord warning label</p> <p>Keyed door locks</p> <p>Magnetic door gasket for positive seal</p> <p>Recessed door handle, freezer; heavy duty door handle, refrigerator</p> <p>Leveling legs</p> <p>Stacking kit</p> <p>LED interior light, Refrigerator</p> <p>Two probe access ports (3/8"); one, refrigerator/one, freezer</p> <p>White powder coated exterior</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Pharmacy refrigerator/freezer toolkit and temperature logs</p> <p>Exterior dimensions: 23 ¾"W x 25 ½"D x 66 ¾"H</p> <p>Two year parts and labor warranty, plus an additional three year compressor parts warranty</p> <p>Manual defrost freezer / cycle defrost refrigerator</p> <p>115V, 60 Hz, 1.74 Amps Refrigerator and 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer</p> <p>UL/C-UL listed</p> <p>Shipping weight: 282 lbs.</p> <p>*Please note this unit requires 2 outlets.</p> |

