

Product Name:	PH-LRP-HC-S16S
Description:	Basic Series Pharmacy/Vaccine Solid Door Refrigerator 16 Cu. Ft. capacity One swing solid door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in sample bottle Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) Calibrated using a NIST traceable device, certificate included Three adjustable shelves (adjustable in ½" increments) Vaccine storage power cord warning label Forced draft air circulation Heavy duty door handle Keyed door lock Magnetic door gasket for positive seal LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (⅜") White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) Pharmacy refrigerator/freezer toolkit and temperature logs One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 25" W x 29 ¾" D x 79" H with casters; 75 ½" H without casters Cycle defrost 115V, 60 Hz, 3 Amps, 1/5 HP ENERGY STAR CERTIFIED ETL listed Shipping weight: 275 lbs.