

Product Name:	PH-LRP-HC-SSP-23
Description:	<p>Plus Series Pharmacy Stainless Steel Refrigerator 23</p> <p>Cu. Ft. Capacity</p> <p>One swing solid door, self-closing, right hinged, not reversible 3 ea.</p> <p>Adjustable Shelves</p> <p>Digital display microprocessor temperature controller Adjustable operating temperature: 36°F to 46°F</p> <p>Calibrated using a NIST traceable device, certificate included</p> <p>Audible and visual Hi/Lo temperature alarm</p> <p>2 probes; 1 in air and 1 in sample bottle</p> <p>Sensor error alarm</p> <p>Continuous product (simulated) temperature display</p> <p>Remote alarm contacts</p> <p>Forced draft air circulation</p> <p>Keyed door lock</p> <p>Magnetic door gaskets for positive seal</p> <p>Integrated door handle</p> <p>LED interior lights are safety shielded and switch controlled Probe access port (3/4")</p> <p>Stainless steel exterior</p> <p>White-powder coated aluminum interior with 304 stainless steel floor Easily cleaned smooth interior</p> <p>Hydrocarbon natural refrigerant (R290)</p> <p>2 year parts and labor warranty, plus an additional 3 year compressor parts warranty</p> <p>4 swivel casters (5 3/4")</p> <p>Exterior dimensions: 26 7/8" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 3 Amps, 1/5 HP</p> <p>ETL listed</p> <p>Shipping weight: 311 lbs.</p>