

Product Name:	PH-LRP-HC-SSP-49G
Description:	Plus Series Pharmacy Stainless Steel Glass Door Refrigerator 49 Cu. Ft. capacity Two swing glass doors, self-closing, not reversible 8 ea. adjustable shelves Digital display microprocessor temperature controller Adjustable operating temperature: 36°F to 46°F) Calibrated using a NIST traceable device, certificate included 2 probes; 1 in air and 1 in sample bottle Sensor error alarm Continuous product (simulated) temperature display Remote alarm contacts Forced draft air circulation Keyed door locks Magnetic door gaskets for positive seal Integrated door handle LED interior lights are safety shielded and switch controlled Probe access port (3/4") Stainless steel exterior White-powder coated aluminum interior with 304 stainless steel floor Easily cleaned smooth interior Hydrocarbon natural refrigerant (R290) 2 year parts and labor warranty, plus an additional 3 year compressor parts warranty 4 swivel casters (5 3/4") Exterior dimensions: 54" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters Cycle defrost 115V, 60 Hz, 4.5 Amps, 1/4 HP ETL listed Shipping weight: 557 lbs.