

Product Name:	PH-LRP-HC-SSP-72G
Description:	<p>Plus Series Pharmacy Stainless-Steel Glass Door Refrigerator</p> <p>72 Cu. Ft. capacity</p> <p>Three swing glass doors, self-closing, not reversible</p> <p>12 ea. adjustable shelves</p> <p>Digital display microprocessor temperature controller</p> <p>Adjustable operating temperature: 36°F to 46°F)</p> <p>Calibrated using a NIST traceable device, certificate included</p> <p>2 probes; 1 in air and 1 in sample bottle</p> <p>Sensor error alarm Continuous product (simulated) temperature display</p> <p>Remote alarm contacts</p> <p>Forced draft air circulation</p> <p>Keyed door locks</p> <p>Magnetic door gaskets for positive seal</p> <p>Integrated door handle</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Probe access port (3/4")</p> <p>Stainless steel exterior</p> <p>White-powder coated aluminum interior with 304 stainless steel floor</p> <p>Easily cleaned smooth interior</p> <p>Hydrocarbon natural refrigerant (R290)</p> <p>2 year parts and labor warranty, plus an additional 3 year compressor parts warranty</p> <p>6 swivel casters (5 3/4")</p> <p>Exterior dimensions: 81" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 7.5 Amps, 1/2 HP</p> <p>ETL listed</p> <p>Shipping weight: 786 lbs.</p>