

Product Name:	PH-LRP-HC-UCBI-0404
Description:	<p> Plus Series Pharmacy/Vaccine Undercounter Refrigerator Built-In  4.6 Cu. Ft. Capacity  One solid swing door, right hinged*  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  2 probes; 1 in air and 1 in sample bottle  Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)  Calibrated using a NIST traceable device, certificate included  Three shelves (two adjustable/one fixed) with guard rail on back to prevent product from being pushed against the back wall  Back wall evaporator cover  Forced draft circulation  Vaccine storage power cord warning label  Recessed door handle  Keyed door lock  Magnetic door gasket for positive seal  Leveling legs  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 23 ¾"W x 24 ½"D x 33 ⅜"H  Cycle defrost  115V, 60 Hz, 1/8 HP, 1.74 Amps  ENERGY STAR CERTIFIED  Energy Consumption 1.15 kWh/day  UL/C-UL listed  Shipping weight: 137 lbs.  *This product is available left hinged, product number PH-LRP-HC-UCBI-0404-LH </p>