

Product Name:	PH-LRP-HC-UCBI-0420-ADA
Description:	Plus Series Pharmacy/Vaccine Undercounter Freezer Built-In 4.2 Cu. Ft. Capacity One solid swing door, right hinged* Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in sample bottle Adjustable operating temperature range: 5°F – -13°F (-15°C – -25°C) Calibrated using a NIST traceable device, certificate included Two fixed shelves Vaccine storage power cord warning label Keyed door lock Magnetic door gasket for positive seal Recessed door handle Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 23 ¾"W x 24 ½"D x 31 15/16"H Manual defrost 115V, 60 Hz, 1.5 Amps, 1/8 HP ADA compliant UL/C-UL listed ENERGY STAR CERTIFIED Shipping weight: 133 lbs. *This product is available left hinged, product number PH-LRP-HC-UCBI-0420-ADA-LH