

Product Name:	PH-LRP-HC-UCFS-0104G-LH
Description:	Plus Series Pharmacy/Vaccine Countertop Refrigerator Freestanding-LEFT HINGED 1 Cu. Ft. Capacity One glass swing door, left hinged Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in glycol bottle Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) Calibrated using a NIST traceable device, certificate included 2 adjustable shelves Back wall evaporator cover Forced draft circulation Vaccine storage power cord warning label Keyed door lock Integrated door handle Magnetic door gasket for positive seal Leveling legs LED interior light Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 17 ¼"W x 19 ¼"D x 21 ¼"H Note: This unit must have 4" clearance on sides and back for adequate ventilation Cycle defrost 115V, 60 Hz, .90 Amps, 1/8 HP ENERGY STAR CERTIFIED Energy Consumption 0.71 kWh/day UL/C-UL listed Shipping weight: 91 lbs.