

Product Name:	LRP-HC-UCBI-0420SS
Description:	<p>Plus Series Undercounter Freezer Built-In  4.2 Cu. Ft. Capacity  One stainless steel swing door, right hinged*  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature range: -15°C – -25°C  Two fixed shelves  Keyed door lock  Recessed door handle  Magnetic door gasket for positive seal  Leveling legs  Probe access port (3/8")  Stainless steel exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 23 ¾" W x 24 ½" D x 33 3/8" H  Manual defrost  115V, 60 Hz, 1.5 Amps, 1/8 HP  ENERGY STAR CERTIFIED  Energy Consumption 0.97 kWh/day  UL/C-UL listed  Shipping weight: 131 lbs.  *This product is available left hinged, product number LRP-HC-UCBI-0420SS-LH</p>