

Product Name:	PH-LRP-HC-UCFS-0430
Description:	<p>Plus Undercounter Pharmacy Freezer Freestanding (-30°C) 4 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: -4°F to -22°F (-20°C to -30°C) Factory Preset at -30°C Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall Bottom footed shelf to allow air circulation Back wall evaporator cover Keyed door lock Magnetic door gasket for positive seal Recessed door handle Leveling legs LED interior light Probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R290), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, compressor parts only Omnicell® and warranty AcuDose RX® compatible Exterior dimensions: 23 3/4" W x 24"D x 32 1/8"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost 115V, 60 Hz, 2.2 Amps, 1/4 HP ENERGY STAR CERTIFIED UL/C-UL listed Shipping weight: 128 lbs.</p>