

Product Name:	LRP-FRP-30
Description:	<p>Plus Series Flammable Storage Refrigerator  30 Cu. Ft. capacity  One solid swing door, right hinged, not reversible  Microprocessor Temperature Controller  Adjustable operating temperature range: 1°-10°C  Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms  Sample simulated digital temperature display; .1° Resolution  Temperature display in Fahrenheit  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Resettable 1-point Min/Max History  Remote alarm contacts  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  Five adjustable shelves  Keyed door lock  Magnetic door gasket for positive seal  Stainless steel door handle  Leveling legs  One probe access port (3/8")  White textured cabinet and door  High density urethane foam cabinet and door insulation  CFC-free refrigerant (R404a), foam insulation and packaging  Two year parts and labor warranty  Exterior dimensions: 35 ¾ "W x 35"D x 78 ½"H  Cycle defrost  115V, 60 Hz, 6 Amps, 1/3 HP  UL/C-UL listed  Units meet NFPA guidelines 45 and 70  Shipping weight: 416 lbs.</p>