



4 Cu. Ft. Capacity Plus Series Undercounter Freezer Freestanding (-20°C)

Model: LRP-HC-UCFS-0520



Description:



The Plus Series Built In Undercounter Freezers feature a microprocessor temperature controller and display, allow for precise temperature control, verification and recovery versus mechanical style thermostats. Critical samples and supplies are security backed with an extensive package of audible and visual alarms.

The Plus Series Free Standing Undercounter Freezers also feature a microprocessor temperature controller and display. The freezer units feature either an automatic or manual defrost operating system. For manual defrost units, periodic defrosting is required.

Freestanding models ("UCFS") require 4" of clearance around back and sides of the unit for proper ventilation. Improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure.

Adjustable operating temperature range from -20°C to -30°C. Factory Preset at -30°C. See "View/Print Product Information" drop down for instructions.

Lab Research Products is Going Green! We are an industry leader in providing a variety of Refrigerators and Freezers that are built with Hydrocarbon Refrigerants, R290 or R600a.

Environmentally friendly, natural hydrocarbon (HC) refrigerants (not containing hydrofluorocarbons) vastly reduce global-warming potential (GWP) while meeting new EPA/SNAP mandates and UL, ASHRAE compliance guidelines.



4 Cu. Ft. Capacity Plus Series Undercounter Freezer Freestanding (-20°C)

Model: LRP-HC-UCFS-0520

Additional Product Images



Details

- 4 Cu. Ft. Capacity
- One solid swing door, right hinged, not reversible
- 2 fixed shelves
- Bottom footed shelf to allow air circulation
- Back wall evaporator cover
- Keyed door lock
- Magnetic door gasket for positive seal
- Recessed door handle
- Leveling legs
- LED interior light
- Probe access port (3/8")
- White cabinet and door
- High density urethane foam cabinet and door insulation
- HFC-free refrigerant (R600a), foam insulation and packaging
- Two year parts and labor warranty, plus an additional three year compressor parts only warranty

Specifications

- Exterior dimensions: 23 5/8" W x 24 1/4"D x 32 3/16"H
- Note: This unit must have 4" clearance on sides and back for adequate ventilation
- Manual defrost
- 115V, 60 Hz, 1.4 Amps, 1/8 HP
- UL/C-UL listed
- Shipping weight: 128 lbs.

Temp Controller

- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- 2 probes; 1 in air and 1 in sample bottle
- Remote alarm contacts
- Adjustable operating temperature range: -15°C – -25°C