



15 Cubic Foot General Purpose Manual Defrost Laboratory Chest Freezers

Model: LRP-MFP-15-C



Description:

The General Purpose Manual Defrost Chest Freezers are also purpose built for critical cold storage applications, but where built-in data logging functionality may not be a necessity.

The all-new temperature display and alarm module offers an array of robust security features with a monitoring and alarm system that is 100% independent of a microprocessor temperature controller.

Periodic manual defrost is required on these units.

Additional Product Images



Details

- 15 Cu. Ft. capacity
- LED Interior light
- Two removable baskets
- Swivel casters (2 3/4")
- Keyed door lock
- Magnetic door gasket for positive seal
- White cabinet and lid
- HFC-free refrigerant (R600a), foam insulation and packaging
- One-year parts and labor warranty, plus an additional four year compressor parts warranty

Specifications

- Exterior dimensions:
 - Width: 56"
 - Depth: 29 1/2"
 - Height: 35 1/8" with casters; 32 3/8" without casters



15 Cubic Foot General Purpose Manual Defrost Laboratory Chest Freezers

Model: LRP-MFP-15-C

(Height with hinged door open 90° is 62 ½")

- Manual defrost
- 115V, 60 Hz, 5 Amps, 1/4 HP
- Shipping weight: 197 lbs.

Temp Controller

- Microprocessor temperature controller:
- Adjustable operating temperature range: -15°C to -25°C
- Standard Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms
- Sample simulated digital temperature display: .1C resolution, convertible to display temperatures in Celsius or Fahrenheit
- Audible and visual high/low temperature alarms
- Sensor error alarm
- Power failure alarm
- Resettable 1-point Min/Max History
- Remote alarm contacts