



23 Cubic Foot Basic Series Pharmacy/Vaccine Solid Door Refrigerator

Model: PH-LRP-HC-S23S



Description:



The Basic Series Solid Door Pharmacy Refrigerators are purpose built for pharmacy and medical applications. The new Basic Series Solid Door Pharmacy Refrigerators feature a digital microprocessor temperature controller allowing for the precise temperature management necessary for critical samples and supplies. A forced air directional refrigeration system provides superior temperature uniformity and recovery after door openings. An automatic cycle-defrost design prevents product damaging temperature spikes. The all-new temperature display and alarm module offers product security with a monitoring and alarm system. Real-time sample simulated temperatures are constantly displayed and supported by audible/visual high and low temperature alarms. Remote alarm dry contacts and an access port are provided to simplify connections to existing facility monitoring systems. Lab Research Products is **Going Green!** We are an industry leader in providing a variety of Refrigerators and Freezers that are built with Hydrocarbon Refrigerants, R290 or R600a. **Environmentally friendly, natural hydrocarbon (HC) refrigerants (not containing hydrofluorocarbons) vastly reduce global-warming potential (GWP) while meeting new EPA/SNAP mandates and UL, ASHRAE compliance guidelines.**

Additional Product Images



Details

- One swing solid door, self-closing, right hinged, not reversible
- Keyed door lock
- Magnetic door gasket for positive seal
- Heavy duty door handle
- Four adjustable shelves (adjustable in ½" increments)
- Vaccine storage power cord warning label
- Forced draft air circulation
- LED interior lights are safety shielded and switch controlled
- Swivel Casters (3 ½")
- Probe access port (?)
- White powder-coated steel exterior and interior



23 Cubic Foot Basic Series Pharmacy/Vaccine Solid Door Refrigerator

Model: PH-LRP-HC-S23S

- Hydrocarbon, natural refrigerant (R290)
- Pharmacy refrigerator/freezer toolkit and temperature logs
- One year parts and labor warranty, plus an additional four year compressor parts warranty
- **Energy Star Certified**

Specifications

- Exterior dimensions: Width: 26 7/8" Depth: 34 3/4" Height: 81 3/4" with casters; 78 1/4" without casters
- Cycle defrost
- 115V, 60 Hz, 3 Amps, 1/5 HP
- ETL Listed
- Shipping weight: 250 lbs

Temp Controller

- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- 2 probes; 1 in air and 1 in sample bottle
- Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)
- Calibrated using a NIST traceable device, certificate included