



Description:



Replaces Part #PH-LRP-S23G The Basic Series Glass Pharmacy Refrigerators are purpose built for pharmacy and medical applications. The new Basic Series Glass Pharmacy Refrigerators feature a digital microprocessor temperature controller allowing for the precise temperature management necessary for critical samples and supplies. A forced air directional refrigeration system provides superior temperature uniformity and recovery after door openings. An automatic cycle-defrost design prevents product damaging temperature spikes. The all-new temperature display and alarm module offers product security with a monitoring and alarm system. Real-time sample simulated temperatures are constantly displayed and supported by audible/visual high and low temperature alarms. Remote alarm dry contacts and an access port are provided to simplify connections to existing facility monitoring systems.

Environmentally friendly, natural hydrocarbon (HC) refrigerants (not containing hydrofluorocarbons) vastly reduce global-warming potential (GWP) while meeting new EPA/SNAP mandates and UL, ASHRAE compliance guidelines

Additional Product Images



Details

- One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible
- Four adjustable shelves (adjustable in ½" increments)
- Vaccine storage power cord warning label
- Forced draft air circulation
- Keyed door lock
- Magnetic door gasket for positive seal
- Heavy duty door handle
- Swivel Casters (3 ½")
- LED interior lights are safety shielded and switch controlled
- Probe access port (?)
- Powder coated steel white interior and exterior
- High density urethane foam cabinet insulation
- Hydrocarbon, natural refrigerant (R290)

- Pharmacy refrigerator/freezer toolkit and temperature logs
- One year parts and labor warranty, plus an additional four year compressor parts warranty
- **Energy Star Certified**

Specifications

- Exterior dimensions: Width: 26 7/8" Depth: 34 3/4" Height: 81 3/4" with casters; 78 1/4" without casters
- Cycle defrost
- 115V, 60 Hz, 3 Amps, 1/5 HP
- ETL listed
- Shipping weight: 336 lbs

Temp Controller

- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- 2 probes; 1 in air and 1 in sample bottle
- Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)
- Calibrated using a NIST traceable device, certificate included