

Product Name:	LRP-FRCB-30
Description:	<p>General Purpose Flammable Storage Refrigerator/Freezer Combination 30 Cu. Ft. Capacity (25 Ref./5 Freezer) One interior door - freezer chamber only/one exterior solid swing door, right-hinged, not reversible Mechanical thermostat Adjustable operating temperature range: 1°C – 10°C Freezer temperature range: -15°C to -20°C Three adjustable refrigerator shelves Keyed door lock Magnetic door gasket for positive seal Stainless steel door handle Leveling legs Probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging One year parts and labor warranty Exterior dimensions: 35 ¾ "W x 35 ½"D x 72 ½ "H Manual defrost 115V, 60Hz, 6 Amps, 1/3 HP UL/C-UL listed Units meet NFPA guidelines 45 and 70 Shipping weight: 396 lbs.</p>