

Product Name:	LRP-FRP-30
Description:	<p>Plus Series Flammable Storage Refrigerator 30 Cu. Ft. capacity One solid swing door, right hinged, not reversible Microprocessor Temperature Controller Adjustable operating temperature range: 1°-10°C Digital Temperature Display and Alarm Module Sample simulated digital temperature display; .1° Resolution Adjustable temperature display in Celsius or Fahrenheit Audible and visual high/low temperature alarms Sensor error alarm Power failure alarm Resettable 1-point Min/Max History Remote alarm contacts Alarm validation (audible, visual, and remote alarm) Alarm mute with ring back Five adjustable shelves Keyed door lock Magnetic door gasket for positive seal Stainless steel door handle Leveling legs One probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging Two year parts and labor warranty Exterior dimensions: 35 ¾ "W x 35"D x 78 ½"H Cycle defrost 115V, 60 Hz, 6 Amps, 1/3 HP UL/C-UL listed Units meet NFPA guidelines 45 and 70 Shipping weight: 416 lbs.</p>