

Product Name:	LRP-HC-49S-TS
Description:	<p>TempLog Plus Series Solid Door Laboratory Refrigerator  49 Cu. Ft. capacity  Two swing solid doors, self-closing, right hinged, not reversible  Microprocessor temperature controller  Adjustable operating temperature range: 1°C - 10°C  Touch Screen Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms  Password-protected control panel  Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit  Active temperature scale bar display  Resettable 1-point min/max history  Digital 24-hour temperature chart  Temperature data logging  USB port for data transfer and software updates  Real-time operational display clock  4-20mA output  Electronic troubleshooting guide  Electronic QR Code to access Owner's Manual &amp; technical support contact information  Adjustable alarm volume  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Door ajar alarm  Low battery alarm  Operation indicator  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Eight adjustable shelves (adjustable in ½" increments)  Forced draft air circulation  Keyed door locks  Magnetic door gasket for positive seal  Heavy duty door handles  LED interior lights are safety shielded and switch controlled  Swivel Casters (3 ½")  Probe access port (¾")  White powder-coated steel exterior and interior  Hydrocarbon, natural refrigerant (R290)  Three year parts and labor warranty, plus an additional two year compressor parts warranty  Exterior dimensions: 54" W x 34 ¾" D x 81 ¾" H with casters; 78 ¼" H without casters  Cycle defrost  115V, 60 Hz, 4.5 Amps, 1/4 HP  ENERGY STAR CERTIFIED  ETL listed  Shipping weight: 446 lbs.</p>